Instruction Manual

Zfx® Cobalt Chromium Eco

Version: 2 / 11. 2024, Replaces version: 1 / 04. 2024

CONTENT

DE (Deutsch) Gebrauchsanweisung Zfx Cobalt Chromium Eco	2
EN (English) Instruction for use Zfx Cobalt Chromium Eco	3
FR (Français) Mode d'emploi Zfx Cobalt Chromium Eco	4
ES (Español) Instrucciones de uso de Zfx Cobalt Chromium Eco	5
IT (Italiano) Istruzioni per l'uso di Zfx Cobalt Chromium Eco	6

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INDIKATION

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Instructions for use

Indications

Beryllium- and Cadmium-free cobalt base blanks are indicated as Type 4 for the manufacture of fixed and removable dentures for crowns and bridge frameworks as well as bridge and implant superstructures according to DIN EN ISO 22674 and EN ISO 9693. The processing is carried out using CAD/CAM technology with suitable milling machines. (Please note manufacturer's specifications).

Contraindications

- In case of known hypersensitivity (allergy) to constituents of the alloy
- all applications not listed under Indications are contraindicated

Technical Specifications:

Type of alloy	4
Density	8,8 g/cm ³
Vickers hardness	285 HV 10/30
CTE values (25-500 °C)	14,4*10 ⁻⁶ K ⁻¹
Oxidising without vacuum	1000 °C
Mass content of all parts .(%)*	Co 63, Cr 29, Mo 6 Mn x, Nb x, Si x, Fe x
Tensile strength	680 MPa
Yield strength	400 MPa
Elongation at break	9 %
E-module	250 GPa

X= Weight-% < 1,0 %

Fabrication

Insert the blank into the machine according to the CAM system operating instructions. and process the disc using tungsten carbide burrs suitable for CoCr. After completion of the milling process, remove the blank with the manufactured frameworks

Separating the frameworks from the blank

The milled restorations can be cut out of the blank using cross-cut tungsten carbide burrs or separating discs suitable for CoCr alloys

Finishing/Cleaning

Machine the separated frameworks with sharp tungsten carbide burrs without high pressure and in one direction only. Avoid sharp edges and interfacing areas. For the desired veneering, the surface must be blasted with aluminium oxid (110 – 180 µm) at a pressure of 3-4 bar, thoroughly evaporated and then degreased with ethanol.

Oxidising

Oxide firing is not necessary for CoCr frameworks. However, if an oxide firing (rapid heating to 1000 °C without vacuum; cool down immediately after reaching the temperature) is carried out to visually check the framework conditioning, the oxide must then be blasted off with disposable aluminium oxide (110 - 180 µm) and then cleaned (e.g. by evaporation).

Veneering with ceramics

All commercially available ceramic materials with a suitable CTE value are suitable for ceramic veneering. The processing instructions recommended by the ceramic manufacturer must be observed, particularly with regard to the cooling conditions after ceramic firing..

Joining (soldering/lasering)

For soldering, a solder and flux suitable for CoCr alloys is recommended. Dimension the soldering block as small as possible and heat it initially at about 600°C for 10 minutes in the oven. Soldering gaps should be between 0,05 and 0,2 mm to avoid solder flowing through solder joints. For laserwelding prior firing, a special laserwelding wire is recommended.

Polishing

After ceramic firing or soldering, any oxide/flux residue must be removed and the framework must be polished to a high gloss using suitable grinding/polishing instruments and polishing paste

Side effects:

- In individual cases, hypersensitivity reactions to constituents of the alloy used cannot be ruled out. The product should not be used in case of known allergies to the constituents of the alloy

Interactions:

- On contact with or in close proximity to alloy of a different type, galvanic effects can occur. In event of sustained local discomfort, other materials must be used for the denture

Safety instructions.

- When designing the prosthetic dentures, ensure adequate wall thickness and connector cross sections.
- These medical devices may only be used and applied by qualified staff in a dental laboratory. Please also observe the material safety data sheet.
- Metal dust is harmful to health. During finishing and stripping, use suitable suction and a respirator mask of type FFP3-EN-2021
- Ensure batch traceability even after the milling process
- The SSCP can be found under the following link: <https://ec.europa.eu/tools/eudamed>
- All serious incidents related to the device shall be reported to the manufacturer and the competent authority of the Member State where the user and/or patient is established

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