CONTENT

EN (English) Instruction Manual GenTek™ Intraoral Scanbody.....	2
BG (Български) Ръководство за употреба Тяло за интраорално сканиране GenTek™.....	4
CS (Česky) Návod k použití Intraorální skenovací tělisko GenTek™.....	6
DA (DANISH) Brugervejledning GenTek™ Intraoral Scan Body.....	8
DE (Deutsch) Gebrauchsanweisung GenTek™ Intraoral Scanbody.....	10
EL (Ελληνικά) Εγχειρίδιο οδηγιών GenTek™ Ενδοστοματικό Scan Body.....	12
ES (Español) Instrucciones de uso GenTek™ Intraoral Scanbody.....	14
ET (Eesti) Kasutusjuhend GenTek™-i suusisene skannimiskeha.....	16
FI (Suomi) Käyttöopas GenTek™ intraoraalinen skannauskappale.....	18
FR (Français) Mode d'emploi Scanbody intra-oral GenTek™.....	20
HR (Hrvatski) Priručnik GenTek™ Nastavak za Intraoralni Skener.....	22
HU (Magyar) Kezelési útmutató GenTek™ szájonbelüli szkennelőtest.....	24
IT (Italiano) Manual d'uso GenTek™ Intraoral Scanbody.....	26
LT (Lietuvių k) Įnstrukcijų vadovas „GenTek™“ skenavimo burnoje kūnas.....	28
LV (Latviski) Lietošanas rokasgrāmata GenTek™ Intraoral Scan Body.....	30
NL (Nederlands) Instructiehandboek Intra-oraal GenTek™ scanlichaam.....	32
NO (Norsk) Instruksjonshåndbok GenTek™ intraoral scanbody.....	34
PL (Polski) Instrukcja obsługi Skaner wewnętrzny GenTek™.....	36
PT (Português) Instruções de uso GenTek™ Intraoral Scanbody.....	38
RO (Română) Manual de instrucțiuni Corp de scanare intraorală GenTek™.....	40
SK (Slovak/Slovakian) Návod na použitie Intraorálne skenovacie teliesko GenTek™.....	42
SL (Slovenščina) Navodila za uporabo Vstavek za intraoralno slikanje GenTek™.....	44
SV (Svenska) Instruktionsmanual GenTek™ Intraoral skanningskropp.....	46

Instruction Manual

GenTek™

Intraoral Scanbody

Version: 7 / 07.2023
Compiled at: 07.07.2023
Effective date: 07.07.2023
Replaces version: 6 / 04.2023

Intended Use

The GenTek™ Intraoral Scanbody is intended to capture the location, angulation and connection orientation of the dental implant in the patient's mouth during the scanning period. The GenTek™ Intraoral Scanbody is available with a genuine connection for the following implant-level Zimmer Biomet Dental implant systems: Certain®, External Hex, TSV®/ Trabecular Metal™, and Eztec®. In addition, abutment level TiBases are available for the Zimmer Biomet Dental Low Profile and Tapered Abutment Systems. For specific dental implant compatibility, please refer to the information on the individual product label. GenTek™ Intraoral Scan Bodies are compatible with Zimmer Biomet Dental prosthetic tooling as outlined in the GenTek™ Components Catalog – please refer to the components catalog located on the Zfx website.

Device description

The GenTek™ Intraoral Scanbody is a two-piece post with a titanium alloy grade 5 ELI (TiAl6V4 ELI) base and a PolyEtherEtherKeton(PEEK) coronal section. The Scanbody engages the connection of the implant and is secured with an integrated screw. The GenTek™ Intraoral Scan Bodies, used in conjunction with the implant-level system, are available in two different heights per package, 4 mm and 7 mm. The GenTek™ Intraoral Scan Bodies used in conjunction with the abutment level system are available in one height per package, the 4 mm.

Product life time

The scan bodies are multiple use with a product life time of 25 autoclave cycles. Product accuracy cannot be guaranteed when exceeding this number of autoclave cycles.

Cleaning and Sterilization

The GenTek™ Intraoral Scanbody is non-sterile. Prior to placement of the GenTek™ Intraoral Scanbody into the patient and prior to each use, the Scanbody must be cleaned and sterilized. Do not disassemble the Scanbody prior to cleaning and sterilization.

Cleaning

1. Rinse the Intraoral Scanbody/Restoration under running cold tap water for 2 minutes.
2. Hand wash instruments in a bath with demineralized water for 5 minutes. Remove debris with a soft brush where possible. Cavities should be rinsed with a syringe. This should be repeated until all debris is removed.
3. Clean the instrument in an ultrasonic bath using a mildly alkaline cleaner (1,5%) at a frequency of at least 40 kHz for 20 minutes.
4. Rinse thoroughly with demineralized water for 1 min.

Sterilization

The recommended sterilization procedure is full cycle pre-vacuum steam sterilization at a temperature of 134°C (273.2°F) for 3 minutes, 4 pulse minimum or full cycle gravity steam sterilization at 134°C(273.2°F) for 3 minutes minimum. Scanbody should not be used prior to 16 minutes after pre-vacuum steam sterilization and 30 minutes after gravity steam sterilization.

Note: Follow the loading instructions issued by the sterilizer manufacturer.

Intra Oral Scanning

- × Remove any provisional abutment, cover screw, healing abutment or similar from the implant upon which the Scanbody is intended to be placed.
- × Confirm that the selected Scanbody is the same size as the implant.
- × Ensure that the seating surface of the implant is clean, free of residue, bone, or soft tissue and is dry when mounting.
- × Affix the GenTek™ Intraoral Scanbody into the implant using a Zimmer Biomet Dental hex driver with a minimum of 5Ncm but no more than 15Ncm maximum torque. In cases with multiple scan bodies, it is important that they are seated with the same torque.
- × Proceed with scanning.

Warnings

- × Do not over-tighten the Scanbody (max torque 15Ncm) as it can cause damage and lead to inaccurate scanning.
- × **DO NOT** acquire the bite scan with the Scanbody still in place. The Scanbody is not designed to withstand biting forces and should be discarded if a patient bites down on it.
- × As the Scanbody is a high precision tool, it may not be used while still warm from the autoclave, as the thermal expansion may lead to inaccurate scanning data.
- × Refer to the instruction for use of the scanner for details on the recommended protocol of the respective scanner.
- × The GenTek™ Intraoral Scanbody may only be used with a compatible digital library file installed in the design software.
- × The GenTek™ Intraoral Scanbody must only be used and handled by dental professionals who have approved training in dentistry.

Instruction Manual

GenTek™

Intraoral Scanbody

Version: 7/07.2023
Compiled at: 07.07.2023
Effective date: 07.07.2023
Replaces version: 6/04.2023

Signs and symbols

Non-sterile

Do not use if package is damaged

Manufacturer Information: Address of legal manufacturer facility

Recommended to fix the Try-In Screw hand-tight

LOT-Number, Symbol followed by Production Date of the lot.

Medical Device

Reference Number, Symbol followed by Manufacturer's Item Number

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Zfx GmbH