



# Zfx™ GenTek™ Pre-milled Abutment Blanks Compatible Abutment Holders

Zfx™ Abutment Holders	Part No.
Zfx™ Pre-milled Abutment Blank Holder, Standard <sup>1</sup>	ZFX02PMSH
Zfx™ Pre-milled Abutment Blank Holder, Multi-Blank for Zfx™ Inhouse5x	ZFX02002069
Zfx™ Pre-milled Abutment Blank Adapter for Glass Ceramic Holder <sup>2</sup>	ZFX02PMA

Medentika® PreFace® Abutment Holders <sup>3</sup>	Part No.
PreFace® Abutment Holder Datron D5	PreFace® Abutment Holders must be ordered directly from the machine manufacturer.
PreFace® Abutment Holder Dental Concept Systems DC1, DC5	
PreFace® Abutment Holder Dental Concept Systems DC5	
PreFace® Abutment Holder imes-icore 140i	
PreFace® Abutment Holder imes-icore 350i, 450i, 550i, 650i, 850i	
PreFace® Abutment Holder imes-icore 750i	
PreFace® Abutment Holder MB Maschinen Cobra Mill 5A1, 5M	
PreFace® Abutment Holder Primacon PFM 24 mediMill	
PreFace® Abutment Holder R+K Organical Multi5X, 5XT-M, Desktop 8	
PreFace® Abutment Holder Röders RXD	
PreFace® Abutment Holder Sirona inLab MC X5	
PreFace® Abutment Holder VHF N4	
PreFace® Abutment Holder VHF S2	
PreFace® Abutment Holder Wissner Gamma 202	

<sup>1</sup> The Zfx™ Pre-milled Abutment Holder, Standard can be used on all machines, which support standard disk holders that are 98 mm in diameter and 10 mm in height.

<sup>2</sup> The Zfx™ Pre-milled Abutment Adapter can be used on all clamping systems that can clamp 10 x 10 mm connection.

<sup>3</sup> Zfx™ GenTek Pre-milled Abutment Blanks can be milled by any milling machine that is equipped with a Medentika® holder. This list is not an all inclusive list of compatible milling machines.